# Dropper® Urine Chemistry Control / Level 1, Level 2

**LOT** 43371  
43372

**REF** 1431-31, 1431-32  
1432-31, 1432-32

2014-09

European Conformity  
CE-Konformitätskennzeichnung  
Conformité aux normes européennes  
Conformità europea  
Conformidad europea

Lot Number  
Bezeichnung  
Designation du lot  
Numero di lotto  
Denominación de lote

Manufactured by  
Hergestellt von  
Fabriqué par  
Fabricado da  
Fabricado por

For in vitro diagnostic use  
In-vitro Diagnosticum  
Pour diagnostic in vitro  
Per uso diagnostico in vitro  
De uso diagnóstico in vitro

Biological Risk  
Biogefährdung  
Risque biologique  
Rischio biologico  
Perigo biológico

Contents of kit  
Inhalt der Packung  
Contenu du coffret  
Contenido della confezione  
Contenido del estuche

Catalog No.  
Bestellnr.  
N° de catalogue  
Catalogo n.  
N° de catálogo

See Product insert  
Siehe Packungsbeilage  
Voir notice d'utilisation  
Vedere il foglietto illustrativo del prodotto  
Consulte el folleto del producto

Authorized Representative  
Bevollmächtigter  
Representant agréé  
Representante autorizzato  
Representante autorizado

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## English

### Intended Use

The Dropper Urine Chemistry Control is intended for use as a control for human urine assay methods. Control materials having known component concentrations are an integral part of diagnostic procedures. Daily monitoring of control values establishes intralaboratory parameters for accuracy and precision of the test method.

### Product Description

The Dropper Urine Chemistry controls are supplied as a ready-to-use liquid requiring no reconstitution. They are prepared from human urine. The Dropper Urine Chemistry controls are fortified to target levels with human proteins and reagent grade chemicals. Preservatives have been added to inhibit microbial growth.

### Warnings and Precautions

**POTENTIAL BIOHAZARDOUS MATERIAL.** All blood donor units comprising the source plasma used in the manufacture of the albumin and globulin used in this product have been tested and found non-reactive for Hepatitis B Surface Antigen, Hepatitis C and HIV 1 & 2 antibody when tested by FDA accepted methods.

No known test method can assure that a product derived from human blood does not contain Hepatitis or HIV virus. It is recommended such samples be handled according to the Centers for Disease Control's Biosafety Level 2 recommendations.

### Storage and Stability

Store the controls at 2° – 8°C. When stored at 2° – 8°C the controls are stable until the expiration date stated on the label. The control has an open vial stability of 18 months from the date of manufacture or until the expiration date. Discard the controls if turbid or if there is any evidence of microbial contamination. Discard controls in the same manner as other biological specimens, according to local guidelines.

### Procedure

Remove the controls from the refrigerator and allow to come to room temperature (18° – 25°C) for 30 to 60 minutes depending on remaining volume. Invert gently to assure homogeneity of the contents. Avoid foaming. Treat the controls as you would a patient sample in accordance with the manufacturer's requirements of the test method. Immediately recap the controls and return to 2° – 8°C when not in use.

### Expected Values

Expected values for the listed lots of controls have been established from interlaboratory data using instrument manufacturers' reagents. Individual laboratory means should fall within the ranges listed. These values should be used as a guide in evaluating the performance of the test methods. Each laboratory should establish its own precision parameters for the methods used to measure each analyte. A quality control program, Quantrol,® is available from Quantimetrix. Please inquire.

Mean values and expected ranges apply to all models of the instrument listed unless otherwise noted.

### Limitations

The expected mean and ranges were established using instrument manufacturer's reagents available at the time of assay. Any future changes made by the manufacturer of a test method may give different values from those previously recovered. Use of methods other than the ones used to establish the expected values may give different values from the ones indicated. Limitations of the test method are included in the package insert for the reagent or instrument being used. Technical updates can be found on our website.

Depending on the instrument and the reagents used to measure creatinine, the mean creatinine values listed may decrease up to 10% over the entire shelf life of the control.

Note: When using Vitros Slides for Creatinine with this control, dilute the control 1:42 instead of the normal 1:21. Run as usual and correct for the dilution.

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