Total Lipids[®] Control

NOTICE: Insert format has changed. Do not discard any pages.
 HINWEIS: Das Format der Packungsbeilage wurde geändert. Bitte keine Seiten wegwerfen.
 AVIS : Le format de la notice d'utilisation a changé. Ne jetez aucune page.
 AVVISO: il formato dell'inserto informativo è stato cambiato. Non gettare via nessuna pagina.
 AVISO: El formato del prospecto ha cambiado. No deseche ninguna página.

Level 1

LOT 74171**REF** Level 1 1741-101**IVD**

2014-04

**ENGLISH**

Intended Use - The Quantimetrix Total Lipids Control is used to monitor the accuracy and precision of HDL Cholesterol, LDL Cholesterol, Total Cholesterol, Triglycerides, Apolipoproteins, and Lp(a) assays.

Product Description - Total Lipids Control is a lyophilized, entirely human serum-based control and is supplied in three levels, 10 x 2 mL vials per level per box. Preservatives have been added to inhibit microbial growth.

HDL Cholesterol/Direct LDL/Total Cholesterol/Triglyceride/Lp(a) values have been assigned using the listed manufacturer's methods on the indicated instruments.

Apolipoprotein values were assigned using the listed manufacturer's immunoturbidimetric methods on the indicated instruments. These methods are based on the International Federation of Clinical Chemistry standardization¹. In the United States, the Apolipoprotein A2, Apolipoprotein C2, Apolipoprotein C3, and Apolipoprotein E methods are for research use only and are not intended for use in diagnostic procedures.

Warnings and Precautions - POTENTIAL BIOHAZARDOUS MATERIAL. All blood donor units comprising the serum pool have been tested and found non-reactive for Hepatitis B surface antigen, Hepatitis C and HIV 1&2 antibody when tested by FDA accepted methods. No known test method can assure that a product derived from human blood does not contain Hepatitis or HIV virus. It is recommended such samples be handled according to the Centers for Disease Control's Bio-Safety Level recommendations.

Storage - When stored at 2-8°C, the unreconstituted material is stable until the expiration date stated on the vial. The reconstituted material is stable for 7 days when stored refrigerated at 2-8°C. Discard the controls if turbid or if there is any evidence of microbial contamination. Discard controls in the same manner as other biological specimens, according to local guidelines.

Procedure - Remove the controls from the refrigerator. Remove stopper and volumetrically add 2.0 mL of purified water. Replace the stopper. Gently swirl to dissolve. Allow the controls to stand 15 minutes. Invert gently and swirl to assure homogeneity of the contents. Let the control stand at least 15 minutes. Swirl gently just prior to use. Avoid foaming. Treat the control as you would a patient sample and test each level in accordance with the reagent manufacturer's requirement of the test method. Recap the control and return it to 2-8°C when not in use.

Expected Values - Values are method dependent. Due to method differences and matrix interactions, the values listed may not always be reproduced by other analytical systems. Each laboratory should establish its own precision parameters. A quality control program, Quantrol[®], is available from Quantimetrix. Please inquire.

Limitations - The degree of matrix sensitivity for many analytical systems is still unknown. Processed materials may exhibit matrix effects which cause them to assay differently from patient specimens on some instruments.² Any changes made by the reagent manufacturer may give different values from the insert range or the laboratory's established ranges. Other limitations of the test method are included in the package insert for the reagent kit or instrument being used.

Note: Occasionally some turbidity may be observed in the controls. This is characteristic of lyophilized controls and is not lipid related. Therefore, the controls, unlike turbid patient samples, should not be diluted. Each control should be run undiluted.

References - ¹Marcovina, S.M., Albers J.J., Standardization of the immunochemical determination of apolipoproteins A1 and B: a report of the IFCC meeting on standardization of apolipoprotein A1 and B measurements. Clin Chem 1989;35:2009-15. ²National Reference System for Cholesterol, U.S. Department of Health and Human Services, Public Health Service, Atlanta GA: Centers for Disease Control.

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Références - ¹Marcovina, S.M., Albers J.J., Standardization of the immunochemical determination of apolipoproteins A1 and B: a report of the IFCC meeting on standardization of apolipoprotein A1 and B measurements. Clin Chem 1989;35:2009-15. ²National Reference System for Cholesterol, U.S. Department of Health and Human Services, Public Health Service, Atlanta GA: Centers for Disease Control.

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Letteratura di riferimento - ¹Marcovina, S.M., Albers J.J., Standardization of the immunochemical determination of apolipoproteins A1 and B: a report of the IFCC meeting on standardization of apolipoprotein A1 and B measurements. Clin Chem 1989;35:2009-15. ²National Reference System

for Cholesterol, U.S. Department of Health and Human Services, Public Health Service, Atlanta GA: Centers for Disease Control.

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Bibliografía - ¹Marcovina, S.M., Albers J.J., Standardization of the immunochemical determination of apolipoproteins A1 and B: a report of the IFCC meeting on standardization of apolipoprotein A1 and B measurements. Clin Chem 1989;35:2009-15. ²National Reference System for Cholesterol, U.S. Department of Health and Human Services, Public Health Service, Atlanta GA: Centers for Disease Control.

GLOSSARY	GLOSSAR	GLOSSAIRE	GLOSSARIO	GLOSARIO
ANALYTES Apolipoprotein Direct HDL Direct LDL Total Cholesterol Triglyceride	ANALYTEN Apolipoprotein Direkt-HDL Direkt-LDL Gesamtcholesterin Triglycerid	ANALYTES Apolipoprotéines HDL direct LDL direct Cholestérol total Triglycérides	ANALITI Apolipoproteina HDL diretto LDL diretto Colesterolo totale Trigliceridi	ANÁLITOS Apolipoproteína HDL directo LDL directo Coolesterol total Triglicéridos
TERMS Instruments Units Level Mean Expected Range	BEGRIFFE Instrumente Einheiten Level Mittelwert Erwarteter Bereich	TERMES Instruments Unités Niveau Moyenne Plage attendue	TERMINI Strumenti Unità Livello Media Range previsto	TÉRMINOS Instrumentos Unidades Concentración Media Intervalo esperado

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English

*Mean and range based on limited data. Each lab should establish its own mean and range.

⁽¹⁾ SI International System of Units

⁽²⁾ Siemens AHDL Flex[®] reagent cartridge, Cat. No. DF48B

⁽³⁾ Kamiya Biomedical Co. K-Assay

⁽⁴⁾ In the United States, the Apolipoprotein A2, Apolipoprotein C2, Apolipoprotein C3, and Apolipoprotein E methods are for research use only and are not intended for use in diagnostic procedures.

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