NOTICE: Insert format has changed. • **HINWEIS:** Das Format der Packungsbeilage wurde geändert. • **AVIS:** Le format de la notice d'utilisation a changé.
AVVISO: il formato dell'inserito informativo è stato cambiato. • **AVISO:** El formato del prospecto ha cambiado.

Total Lipids® Control / Level 1

LOT 74181 **REF** 1741-101 **2016-02** **CE** **IVD**

CE European Conformity CE-Konformitätskennzeichnung Conformité aux normes européennes Conformità europea Conformidad europea	LOT Lot Number Bezeichnung Designation du lot Numero di lotto Denominación de lote	IVD Manufactured by Hergestellt von Fabriqué par Fabricato da Fabricado por	IVD For in vitro diagnostic use In-vitro Diagnostikum Pour diagnostic in vitro Per uso diagnostico in vitro De uso diagnóstico in vitro	Cont. Biological Risk Biogefährdung Risque biologique Rischio biologico Peligro biológico	Cont. Contents of kit Inhalt der Packung Contenu du coffret Contenuto della confezione Contenido del estuche	REF Catalog No. Bestellnr. N° de catalogue Catalogo n. N° de catálogo	! See Product Insert Siehe Packungsbeilage Voir notice d'utilisation Vedere il foglietto illustrativo del prodotto Consulte el folleto del producto	EC REP Authorized Representative Bevollmächtigter Représentant agréé Rappresentante autorizzato Representante autorizado	! Temperature limitation Temperaturbegrenzungen Limites de température limiti di temperatura límite de temperatura	! Use by (last day of month) Verwendbar bis (letzter Tag des Monats) Utilisable jusqu'à (dernier jour du mois indiqué) Da utilizzare prima del (ultimo giorno del mese) Estable hasta (último día del mes)
--	--	---	---	---	--	---	---	--	--	--

English

Intended Use

The Quantimetrix Total Lipids Control is used to monitor the accuracy and precision of HDL Cholesterol, LDL Cholesterol, Total Cholesterol, Triglycerides, Apolipoproteins, and Lp(a) assays.

Product Description

Total Lipids Control is a lyophilized, entirely human serum-based control and is supplied in three levels, 10 x 2 mL vials per level per box. Preservatives have been added to inhibit microbial growth.

HDL Cholesterol/Direct LDL/Total Cholesterol/Triglyceride/Lp(a) values have been assigned using the listed manufacturer's methods on the indicated instruments.

Apolipoprotein values were assigned using the listed manufacturer's immunoturbidimetric methods on the indicated instruments. These methods are based on the International Federation of Clinical Chemistry standardization¹. In the United States, the Apolipoprotein A2, Apolipoprotein C2, Apolipoprotein C3, and Apolipoprotein E methods are for research use only and are not intended for use in diagnostic procedures.

Warnings and Precautions

POTENTIAL BIOHAZARDOUS MATERIAL. All blood donor units comprising the serum pool have been tested and found non-reactive for Hepatitis B surface antigen, Hepatitis C and HIV 1&2 antibody when tested by FDA accepted methods. No known test method can assure that a product derived from human blood does not contain Hepatitis or HIV virus. It is recommended such samples be handled according to the Centers for Disease Control's Bio-Safety Level recommendations.

Storage

When stored at 2-8°C, the unreconstituted material is stable until the expiration date stated on the vial. The reconstituted material is stable for 7 days when stored refrigerated at 2-8°C. Discard the controls if turbid or if there is any evidence of microbial contamination. Discard controls in the same manner as other biological specimens, according to local guidelines.

Procedure

Remove the controls from the refrigerator. Remove stopper and volumetrically add 2.0 mL of purified water. Replace the stopper. Gently swirl to dissolve. Allow the controls to stand 15 minutes. Invert gently and swirl to assure homogeneity of the contents. Let the control stand at least 15 minutes. Swirl gently just prior to use. Avoid foaming. Treat the control as you would a patient sample and test each level in accordance with the reagent manufacturer's requirement of the test method. Recap the control and return it to 2-8°C when not in use.

Expected Values

Values are method dependent. Due to method differences and matrix interactions, the values listed may not always be reproduced by other analytical systems. Each laboratory should establish its own precision parameters. A quality control program, Quantrol®, is available from Quantimetrix. Please inquire.

Limitations

The degree of matrix sensitivity for many analytical systems is still unknown. Processed materials may exhibit matrix effects which cause them to assay differently from patient specimens on some instruments.² Any changes made by the reagent manufacturer may give different values from the insert range or the laboratory's established ranges. Other limitations of the test method are included in the package insert for the reagent kit or instrument being used.

Note: Occasionally some turbidity may be observed in the controls. This is characteristic of lyophilized controls and is not lipid related. Therefore, the controls, unlike turbid patient samples, should not be diluted. Each control should be run undiluted.

References

¹Marcovina, S.M., Albers J.J., Standardization of the immunochemical determination of apolipoproteins A1 and B: a report of the IFCC meeting on standardization of apolipoprotein A1 and B measurements. Clin Chem 1989;35:2009-15. ²National Reference System for Cholesterol, U.S. Department of Health and Human Services, Public Health Service, Atlanta GA: Centers for Disease Control.

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English

*Mean and range based on limited data. Each lab should establish its own mean and range.
 (1) SI International System of Units
 (2) Siemens AHDL Flex® reagent cartridge, Cat. No. DF48B
 (3) Kamiya Biomedical Co. K-Assay
 (4) In the United States, the Apolipoprotein A2, Apolipoprotein C2, Apolipoprotein C3, and Apolipoprotein E methods are for research use only and are not intended for use in diagnostic procedures.

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