Lipisure® Serum Lipoprotein Control / Level 1

LOT 95060 **REF** 48-7060 **2014-10** **CE** **IVD**

CE European Conformity CE-Konformitätskennzeichnung Conformité aux normes européennes Conformità europea Επιστημονική συσκευασία
LOT Lot Number Bezeichnung Désignation du lot Numero di lotto Determinación de lote Αρ. αριθμολογίας
IVD Manufactured by Ηergestellt von Fabriqué par Fabricato da Fabricado por Παρασκευαστής από την
Cont. For in vitro diagnostic use In-vitro Diagnostikum Pour diagnostic in vitro Per uso diagnostico in vitro De uso diagnostico in vitro Για In Vitro (συσκευασία) χρήση
REF Biological Risk Biologische Gefahr Risiko biologique Ρίσκο biologico Περιγινωτικό κίνδυνος
Cont. Contents of kit Inhalt der Packung Contenu du coffret Contenido della confezione Contenido del estuche Περιεχόμενα του σετ
REF Catalog No. Bestellnr. N° de catalogue Catalogo n. N° de catálogo №, артикула
See Product Insert Siehe Packungsbeilage Voir notice d'utilisation Vedere il foglietto illustrativo del prodotto Consulte el folleto del producto Βλ. οδηγό προϊόντος
EC REP Authorized Representative Bevollmächtigter Représentante agréé Representante autorizado Representante autorizado Εξουσιοδοτημένος εκπαιδευτικός
Temperature limitation Temperaturbegrenzungen Limites de température limiti di temperatura limite de temperatura Ορίσματα θερμοκρασίας
Use by (last day of month) Verwendbar bis (letzter Tag des Monats) Utilisable jusqu'à (dernier jour du mois indiqué) Da utilizzare prima del (ultimo giorno del mese) Estable hasta (ultimo día del mes) Αρτέψ (εγκαιρίως) μέχρι του μήνα

English

Intended Use and Indications for Use

The Quantimetrix Lipisure Lipoprotein Control is intended for monitoring lipoprotein fractions (VLDL to HDL) as well as subfractions (Mid-bands, LDL subfractions 1, 2 and 3) using the Lipoprint LDL System to validate the separation and quantitation of lipoproteins in patient serum samples. Control materials with known component concentrations are an integral part of diagnostic procedures. Daily monitoring of control values establishes inter- and intralaboratory parameters for accuracy and precision of the test method.

Product Description

Lipisure is prepared from human serum and supplied lyophilized. It is assayed for VLDL, Midbands C, B and A, LDL subfractions, HDL and Total LDL using the Quantimetrix Lipoprint LDL System. The lyophilized material contains stabilizers and preservatives. The control is provided as a set of 4 x 0.5 ml vials (0.5 ml after reconstitution).

FOR IN VITRO DIAGNOSTIC USE ONLY

Warnings and Precautions

POTENTIAL BIOHAZARDOUS MATERIAL. The Lipisure Control is prepared in a human serum matrix and contains other human source materials. All blood donor units comprising the serum pool have been tested and found nonreactive for Hepatitis B surface antigen, Hepatitis C and HIV 1 & 2 antibodies when tested by FDA accepted methods.

No known test method can assure that a product derived from human source materials does not contain Hepatitis or HIV virus. Handle these controls with the same precautions used when handling any potentially infectious material. It is suggested that such samples be handled according to the Centers for Disease Control's Bio-Safety Level 2 recommendations.

Storage and Stability

- When stored at 2-8° C the unopened, lyophilized material is stable until the expiration date stated on the vial label.
- The reconstituted material is stable for 5 days when kept at 2-8° C.
- Discard the control if it is turbid or if it exhibits any signs of microbial contamination.

Materials Needed but Not Provided

- Deionized or distilled water
- A 0.5 mL pipet

Procedure

- Reconstitute the control with 0.5 mL of purified water. Gently swirl to dissolve.
- Allow the reconstituted control to stand for at least 10 minutes at ambient temperature before use.
- Invert gently to assure complete homogeneity of the solution. Avoid foaming.
- Handle the control the same way you would a patient sample. Use as indicated in the Lipoprint LDL System package insert.
- After use immediately recap the control and return it to refrigerated storage (2-8° C).
- Include at least one control sample with each run of patient samples.

Quantitative Results

The control will generate a normal lipoprotein profile when used with the Lipoprint LDL System. It will typically exhibit bands for VLDL, Midbands-C, B and A (these comprise IDL), LDL-1, 2 and 3 and HDL when using the Lipoprint LDL template. See Figure 1 below.


Figure 1. Example of Lipoprint band pattern for Lipisure Control using the template

Quantitative Results

The gel tubes are scanned and the relative area of each lipoprotein subfraction is determined from the densitometric scan as described in the Lipoprint LDL System package insert. A representative scan is depicted in Figure 2.


Figure 2. Example of a typical control profile

Expected Values

The assigned values have been established from repetitive runs of the control on the Lipoprint LDL System. See Table on reverse side of insert.

Limitations

- Use of this control with other electrophoresis test systems or precipitation methods for the determination of LDL and HDL may give values other than those assigned for the Lipoprint LDL System.
- Lipisure is not intended for use as a control for total cholesterol or triglycerides.
- Lipisure Level 1 is not intended for use as a control for LDL subfractions 5-7.

Assigned values by Lipoprint LDL System

	Units mg/dL ²		SI Units ¹ mmol/L ³	
	Mean	Range	Mean	Range
VLDL	17	11 – 22	0.440	0.284 – 0.569
Mid-C	12	8 – 15	0.310	0.207 – 0.388
Mid-B	7	5 – 8	0.181	0.129 – 0.207
Mid-A	6	4 – 8	0.155	0.103 – 0.207
LDL-1	19	17 – 22	0.491	0.440 – 0.569
LDL-2	14	10 – 17	0.362	0.259 – 0.440
LDL-3	2	0 – 5	0.052	0.000 – 0.129
LDL-4	0	0 – 1	0.000	0.000 – 0.026
HDL	23	20 – 26	0.595	0.517 – 0.672
Total LDL ⁴	60	55 – 65	1.552	1.422 – 1.681

¹SI International System of Units
²Based on total cholesterol of 100 mg/dL
³Based on total cholesterol of 2.59 mmol/L
⁴Total LDL equals the sum of all Midbands and LDL subfractions present

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