**NOTICE:** Insert format has changed. • **HINWEIS:** Das Format der Packungsbeilage wurde geändert. • **AVIS:** Le format de la notice d'utilisation a changé.  
**AVVISO:** il formato dell'insero informativo è stato cambiato. • **AVISO:** El formato del prospecto ha cambiado.

# Ammonia Alcohol Control / Level 1, 2 & 3

**LOT** 31391  
31392  
31393

**REF** 1311-31  
1312-31  
1313-31

2014-01

European Conformity  
CE-Konformitätskennzeichnung  
Conformité aux normes européennes  
Conformità europea  
Conformidad europea

Lot Number  
Bezeichnung  
Désignation du lot  
Numero di lotto  
Denominación de lote

Manufactured by  
Hergestellt von  
Fabriqué par  
Fabbricato da  
Fabricado por

For in vitro diagnostic use  
In-vitro Diagnosticum  
Pour diagnostico in vitro  
Per uso diagnostico in vitro  
De uso diagnostico in vitro

Biological Risk  
Biosigefährdung  
Risque biologique  
Rischio biologico  
Peligro biológico

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## English

### Intended Use

The Quantimetrix Ammonia/Alcohol control is intended as a means of monitoring various ammonia and ethanol assay methods to validate quantitation of patient samples.

Control materials having known component concentrations are an integral part of diagnostic procedures. Daily monitoring of control values establishes intralaboratory parameters for accuracy and precision of the test method.

### Product Description

The Quantimetrix Ammonia/Alcohol Control is supplied in three levels, 3 x 5 mL, as a ready-to-use liquid requiring no reconstitution or dilution. It is fortified with ethanol and reagent grade chemicals. Preservatives have been added to inhibit microbial growth.

### Precautions

For in vitro diagnostic use only. Observe the same safety precautions used when handling potentially infectious biological material. Dispose of properly.

### Storage and Stability

Store the control at 2 – 8° C. When stored unopened at 2 – 8° C, the control is stable until the expiration date stated on the label. The control has an open vial stability of four months after opening and use. Discard the control if it is turbid or there is any evidence of microbial contamination. Discard controls in the same manner as other biological specimens, according to local guidelines.

### Procedure

Remove the controls from the refrigerator and allow to come to room temperature (18 – 25° C), at least 15 minutes, depending on remaining volume. Gently invert the control to assure homogeneity of the contents. Avoid foaming. Treat the control as you would a patient sample in accordance with the manufacturer's requirements of the test method. Immediately recap the control and return it to 2 – 8° C when not in use.

### Expected Range

The expected ranges have been established from interlaboratory data. Each laboratory should establish its own precision parameters. A quality control program, Quantrol,™ is available from Quantimetrix. Please inquire.

### Limitations

Any future changes made by the manufacturer of a test method may give different values from the indicated range. Limitations of the test method are included in the package insert for the reagent kit or instrument being used. Technical updates can be found on our website.

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