NOTICE: Insert format has changed. • **HINWEIS:** Das Format der Packungsbeilage wurde geändert. • **AVIS:** Le format de la notice d'utilisation a changé.
AVVISO: il formato dell'insero informativo è stato cambiato. • **AVISO:** El formato del prospecto ha cambiado.

Bilirubin Control – Pediatric / Level 1 & 2

LOT 33421
33422

REF 1331-31
1332-31

2013-12

European Conformity
CE-Konformitätskennzeichnung
Conformité aux normes européennes
Conformità europea
Conformidad europea

Lot Number
Bezeichnung
Désignation du lot
Número de lotto
Denominación de lote

Manufactured by
Hergestellt von
Fabriqué par
Fabricato da
Fabricado por

In vitro diagnostic use
In-vitro Diagnosticum
Pour diagnostic in vitro
Per uso diagnostico in vitro
De uso diagnóstico in vitro

Biological Risk
Biosgefahr
Risque biologique
Rischio biologico
Peligro biológico

Contents of kit
Inhalt der Packung
Contenu du coffret
Contenido della confezione
Contenido del estuche

Catalog No.
Besteller
N° de catalogue
Catalogo n.
N° de catálogo

See Product Insert
Siehe Packungsbeilage
Voir notice d'utilisation
Vedere il foglietto illustrativo del prodotto
Consulte el folleto del producto

Authorized Representative
Bevollmächtigter
Représentant agréé
Representante autorizzato
Representante autorizado

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English

Intended Use

The Quantimetrix Pediatric Bilirubin Control is intended as a means of monitoring various bilirubin assay methods to validate quantitation of patient samples.

Control materials having known component concentrations are an integral part of diagnostic procedures. Daily monitoring of control values establishes intralaboratory parameters for accuracy and precision of the test method.

Product Description

The Quantimetrix Pediatric Bilirubin Control is supplied as a liquid, 3 x 3 mL vials per box. It is prepared from purified bilirubin in a human protein base. Stabilizers and preservatives have been added.

Warnings and Precautions

POTENTIAL BIOHAZARDOUS MATERIAL. All blood donor units comprising the protein pool have been tested and found non-reactive for Hepatitis B surface antigen, and for HCV antibody and HIV antibody when tested by FDA accepted methods. No known test method can assure that a product derived from human blood does not contain Hepatitis or HIV virus. It is suggested that such samples be handled according to the Centers for Disease Control's Bio-Safety Level 2 recommendations.

Storage and Stability

The control should be stored at 2 – 8° C. When stored unopened at 2 – 8° C, the control is stable until the expiration date stated on the label. If the rubber stopper is removed or air is introduced into the vial, the stability of the control is limited to five (5) days at 2 – 8° C. **Longer stability may be achieved by removing the sample with a syringe.** Discard the control if it is turbid, if there is any evidence of microbial contamination or if it exhibits any green color. Discard controls in the same manner as other biological specimens, according to local guidelines.

Procedure

Remove the control from the refrigerator and allow it to come to room temperature (18 – 25° C) for 10 to 15 minutes. Gently invert the control to assure homogeneity of the contents. Avoid foaming the control. For maximum stability, lift up the center tab of the metal seal without removing the seal. Pierce the rubber stopper with a syringe and withdraw the required quantity of control. **Do not introduce air into the vial.** Do not return unused material to the vial. Treat the control as you would a patient sample in accordance with the manufacturer's requirements of the test method.

Expected Range

The expected ranges have been established from interlaboratory data. Each laboratory should establish its own precision parameters. A quality control program, Quantrol™, is available from Quantimetrix. Please inquire.

Limitations

Bilirubin is sensitive to light and oxygen. Exposure to these factors may decrease the values from those stated in the expected ranges.

Any future changes made by the manufacturer of a test method may give different values from the indicated range. Limitations of the test method are included in the package insert for the reagent kit or instrument being used. Technical updates can be found on our website.

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Beckman Coulter, Inc., Synchron®
Beckman Coulter, Inc., AU400®, AU640® and AU2700®
Roche Diagnostics Corporation, Hitachi, Cobas Integra®, Cobas Mira®
Advanced Instruments, Inc., Bilirubin Stat-Analyzer Model BR2