Complete D 25-OH Vitamin D Control / Level 1, Level 2

LOT 29081
29082

REF 1291-31
1292-31

2014-10

European Conformity
CE-Konformitätskennzeichnung
Conformité aux normes européennes
Conformità europea
Conformidad europea

Lot Number
Bezeichnung
Designation du lot
Numero di lotto
Denominación de lote

Manufactured by
Hergestellt von
Fabrique par
Fabricato da
Fabricado por

For in vitro diagnostic use
In-vitro Diagnosticum
Pour diagnostic in vitro
Per uso diagnostico in vitro
De uso diagnóstico in vitro

Biological Risk
BIOGefährdung
Risque biologique
Rischio biologico
Peligro biológico

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English

Intended Use

The Complete D 25-OH Vitamin D control is intended to monitor the performance of clinical assays used to quantitate total 25-OH Vitamin D.

Use of quality control materials is an integral part of diagnostic procedures. Daily monitoring of control values establishes intralaboratory parameters for accuracy and precision of the test method.

Product Description

The Complete D 25-OH Vitamin D control is supplied in two levels, 3 mL per vial. Each kit contains 3 vials of level 1 or level 2 controls. The controls are liquid, ready to use and require no reconstitution or dilution (See limitations). The controls are prepared in a human serum matrix fortified with analytes and preservatives. This control contains both 25-OH Vitamin D2 and D3; and when combined represent clinically relevant levels of Total Vitamin D.

Warnings and Precautions

POTENTIAL BIOHAZARDOUS MATERIAL. The Complete D Controls are prepared in a human serum matrix and contain other human source materials. All blood donor units comprising the serum pool have been tested and found nonreactive for Hepatitis B surface antigen, Hepatitis C and HIV 1 & 2 antibodies when tested by FDA accepted methods.

No known test method can assure that a product derived from human source materials does not contain Hepatitis or HIV virus. Handle these controls with the same precautions used when handling any potentially infectious material. It is suggested that such samples be handled according to the Centers for Disease Control's Bio-Safety Level 2 recommendations.

Storage and Stability

Store the control at 2 – 8°C. When stored unopened at 2 – 8°C, the control is stable until the expiration date stated on the label. The control has an open vial stability of one year when stored at 2 – 8°C and 30 days when stored at room temperature (18 – 25°C).

Procedure

Remove the controls from the refrigerator and allow to come to room temperature (18 – 25°C), at least 15 minutes, depending on remaining volume. Gently invert the control to assure homogeneity of the contents. Avoid foaming. Treat the control as you would a patient sample in accordance with the manufacturer's requirements of the test method. Immediately recap the control when not in use.

Expected Range

The expected ranges have been established from interlaboratory data. Each laboratory should establish its own precision parameters. A quality control program, Quantrol,® is available from Quantimetrix. Please inquire.

Limitations

Note: For Diasorin users, dilute Level 1 and Level 2 by 75% (3 parts control to 1 part Liaison® Vitamin D Specimen Diluent). The expected mean and ranges were established using instrument manufacturer's reagents available at the time of assay. Any future changes made by the manufacturer of a test method may give different values from those previously recovered. Use of methods other than the ones used to establish the expected values may give different values from the ones indicated. Limitations of the test method are included in the package insert for the reagent or instrument being used. Technical updates can be found on our website.

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