Microbumin® Microalbumin Control / Level 1 & 2

LOT 25331
25332 **REF** 1250-01 2014-05 **IVD**

European Conformity
CE-Konformitätskennzeichnung
Conformità ai norme europee
Conformità europea
Conformidad europea

Lot Number
Bezeichnung
Designation du lot
Numero di lotto
Denominación de lote

Manufactured by
Hergestellt von
Fabriqué par
Fabricato da
Fabricado por

For in vitro diagnostic use
In-vitro Diagnosticum
Pour diagnostic in vitro
Per uso diagnostico in vitro
De uso diagnostico in vitro

Biological Risk
Biogefährdung
Risque biologique
Rischio biologico
Peligro biológico

Contents of kit
Inhalt der Packung
Contenu du coffret
Contenido della confezione
Contenido del estuche

Catalog No.
Bestellnr.
N° de catalogue
Catalogo n.
N° de catálogo

See Product Insert
Siehe Packungsbeilage
Voir notice d'utilisation
Vedere il foglietto illustrativo del prodotto
Consulte el folleto del producto

Authorized Representative
Bevollmächtigter
Représentante agréé
Representante autorizzato
Representante autorizado

Temperature Limitation
Temperaturbegrenzungen
Limites de température
limiti di temperatura
límite de temperatura

Use by (last day of month)
Verwendbar bis (letzter Tag des Monats)
Utilisable jusqu'à (dernier jour du mois indiqué)
Da utilizzare prima del (ultimo giorno del mese)
Estable hasta (ultimo día del mes)

English

Intended Use

The Quantimetrix Urine Microalbumin Controls are intended as a means of monitoring various microalbumin assay methods to validate quantitation of patient samples. Control materials having known component concentrations are an integral part of diagnostic procedures. Daily monitoring of control values establishes intralaboratory parameters for accuracy and precision of the test method.

Product Description

The Quantimetrix Urine Microalbumin Controls are supplied in two levels, 5 x 7 mL each level per box. They are ready to use, liquid, requiring no reconstitution. They are prepared from human urine, fortified to target levels with human albumin and creatinine. Preservatives including sodium azide have been added to inhibit microbial growth.

Warnings and Precautions

POTENTIAL BIOHAZARDOUS MATERIAL. All blood donor units comprising the source plasma used in the manufacture of the albumin have been tested and found non-reactive for Hepatitis B Surface Antigen and HIV antibody when tested by FDA accepted methods. No known test method can assure that a product derived from human blood does not contain Hepatitis or HIV virus. It is recommended that such samples be handled according to the Center for Disease Control's Bio-Safety Level 2 recommendations. Dispose of carefully. Sodium azide may accumulate in plumbing traps and pose a threat of explosion.

Storage and Stability

The controls should be stored at 2-8° C. When stored at 2-8° C, the controls are stable until the expiration date stated on the label. When stored at 2-8° C between each use, the controls are stable for six months after opening. When using the control with the Micral Test Strips, the control will remain stable for six months after opening, until after 10 uses, or until the expiration date, whichever occurs first. Discard the controls if turbid or if there is any evidence of microbial contamination. Discard controls in the same manner as other biological specimens, according to local guidelines.

Procedure

Remove the controls from the refrigerator and allow to come to room temperature (20-25° C), about 15 - 30 minutes. Invert gently to assure homogeneity of the contents. Avoid foaming. Treat the control as you would a patient sample in accordance with the manufacturer's requirements of the test method. Immediately recap the controls and return to 2-8° C when not in use.

Expected Values

Expected values for the listed lots of controls have been established from interlaboratory data using instrument manufacturers' reagents. Individual laboratory means should fall within the ranges listed. These values should be used as a guide in evaluating the performance of the test methods. Each laboratory should establish its own precision parameters for the methods used to measure each analyte. A quality control program, Quantrol,™ is available from Quantimetrix. Please inquire.

Mean values and expected ranges apply to all models of the instrument listed unless otherwise noted.

Limitations

The expected mean and ranges were established using instrument manufacturer's reagents available at the time of assay. Any future changes made by the manufacturer of a test method may give different values from those previously recovered. Use of methods other than the ones used to establish the expected values may give different values from the ones indicated. Limitations of the test method are included in the package insert for the reagent or instrument being used.

Depending on the instrument and the reagents used to measure creatinine, the mean creatinine values listed may decrease up to 10% over the entire shelf life of the control. NOTE: For Creatinine assays using the Dimension and Vitros slides, dilute the Level 2 control with an equal volume of diluent (1 part control to 1 part enzyme diluent). Run as usual and correct for the dilution.

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