Technical update: Spinalscopics® Spinal Fluid Cell Count Control

Tech number: 320-12-46-1 • Users: All Spinalscopics Users • Date issued: 11/16/12

English

To provide greater clarity and reduce variability, instructions for the Spinalscopics product has been updated as follows:

Revised Hemocytometer Cell Count Procedure

1. Remove the control from the refrigerator and replace the cap on the control bottle with the bulb dropper cap included in the control box.
2. Allow the control to come to room temperature (**18-25°C**) for at least 15 minutes.
3. Mix the control thoroughly by inverting the bottle several times and by squeezing the bulb in the cap, and by aspirating and expelling the control through the glass dropper attached to the cap (**at least 10 times**) **immediately** prior to use to ensure homogeneity of the content. Thorough mixing with each use is important in order to obtain reproducible results. Avoid foaming.

Note: To prevent cells from settling and to ensure homogeneity of sample, the hemocytometer chambers must be charged immediately after thorough mixing. Do not allow cells to settle in bottle prior to charging hemocytometer chambers.

4. Using the glass dropper provided, **charge both sides of the hemocytometer chamber by placing the glass dropper tip alongside cover-slip for each chamber. To prevent overflow chambers, gently squeeze the bulb to slowly charge the chamber until full.**
5. Immediately recap the controls after each use. The controls can be stored for one month at room temperature when not in use. Treat the controls as you would a patient sample.
6. Allow the cells to settle by placing the charged hemocytometer in a humidified chamber for 10 minutes before counting.
7. Count the cells in all nine (9) large squares on both sides of the hemocytometer. **If there are cells on the line, count only those on the top and right side of each square.**
8. Average the number of cells counted on both sides of chamber. Calculate the total number of erythrocytes (RBC) and lymphocytes (WBC) per μL as follows:

Average cell count (cells/ μL) \div 0.9 = Total cell count (cells/ μL)

Note: Clinical data demonstrates that the majority of leukocytes (WBC) found in Cerebrospinal Spinal Fluid (CSF) are lymphocytes. To more closely match clinical sample, the leukocytes (WBC) population in this product has been lymphocytes enriched. As lymphocytes are similar in size to erythrocytes (RBC), it is important to ensure that laboratory personnel can differentiate lymphocytes (WBC) from erythrocytes (RBC).

Any deviation from this procedure stated above may introduce variables that can result in values outside specified product insert range. Each laboratory should establish its own precision parameters.

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