Sweat Control

for Cystic Fibrosis Testing

Level **1,2&3**
CE

LOT 19191
19192
19193

REF 1190-01

IVD

2014-08

| Analytes / Instruments | Units | Level 1 - 19191 | | Level 2 - 19192 | | Level 3 - 19193 | |
|---|--------|-----------------|----------------|-----------------|----------------|-----------------|----------------|
| | | Mean | Expected Range | Mean | Expected Range | Mean | Expected Range |
| Chloride | | | | | | | |
| Chloridometer | mmol/L | 25 | 20 - 30 | 49 | 39 - 59 | 93 | 75 - 112 |
| ISE Electrode | | | | | | | |
| Beckman Synchron | mmol/L | 26 | 21 - 32 | 51 | 41 - 61 | 96 | 77 - 116 |
| Dade Dimension | mmol/L | 30 | 24 - 36 | 55 | 44 - 66 | 103 | 82 - 123 |
| Ortho Vitros* | mmol/L | n/a | <50 | 51 | <50 - 61 | 96 | 77 - 115 |
| Roche Hitachi* | mEq/L | n/a | <80 | n/a | <80 | 91 | <80 - 109 |
| Potassium | | | | | | | |
| ISE Electrode | | | | | | | |
| Beckman Synchron | mmol/L | 9 | 7 - 10 | 18 | 14 - 21 | 34 | 27 - 41 |
| Dade Dimension | mmol/L | 8 | 7 - 10 | 17 | 14 - 21 | 34 | 27 - 40 |
| Ortho Vitros* | mmol/L | 9 | 7 - 11 | 18 | 15 - 22 | 37 | 29 - 44 |
| Roche Hitachi* | mEq/L | 8 | 6 - 10 | 18 | 14 - 21 | 37 | 29 - 44 |
| Sodium | | | | | | | |
| ISE Electrode | | | | | | | |
| Beckman Synchron | mmol/L | 27 | 22 - 33 | 54 | 43 - 65 | 103 | 82 - 124 |
| Dade Dimension | mmol/L | 29 | 23 - 35 | 54 | 44 - 65 | 100 | 80 - 121 |
| Ortho Vitros* | mmol/L | 24 | 19 - 29 | 51 | 41 - 61 | 105 | 84 - 126 |
| Roche Hitachi* | mEq/L | 32 | 26 - 38 | 56 | 44 - 67 | 104 | 83 - 125 |
| Conductivity | | | | | | | |
| Wescor Sweat-Chek™* | mmol/L | 37 | 30 - 45 | 70 | 56 - 84 | 124 | 99 - 149 |
| Osmolality | | | | | | | |
| Advanced® Series Osmometers | mOs/kg | 73 | 58 - 87 | 143 | 114 - 171 | 272 | 217 - 326 |
| VAPRO® Vapor Pressure Osmometer* | mOs/kg | 68 | 54 - 82 | 139 | 111 - 167 | 265 | 212 - 318 |
| N/A Not applicable *Mean and range based on limited data. | | | | | | | |
| Advanced Instruments, Inc.; Beckman Coulter Inc.; Siemens, Inc.; Instrumentation Laboratories; Ortho-Clinical Diagnostics; Wescor, Inc. | | | | | | | |

| GLOSSARY | GLOSSAR | GLOSSAIRE | GLOSSARIO | GLOSARIO |
|---|--|---|---|--|
| ANALYTES Chloride Potassium Sodium Conductivity Osmolality | ANALYTEN Chlor Kalium Natrium Leitfähigkeit Osmolalität | ANALYTES Chlorure Potassium Sodium Conductivité Osmolalité | ANALITI Cloruro Potassio Sodio Conduktività Osmolalità | ANÁLITOS Cloruro Potasio Sodio Conductividad Osmolalidad |
| TERMS Instruments Units Level Mean Expected Range Chloridometer ISE Electrode | BEGRIFFE Instrumente Einheiten Level Mittelwert Erwarteter Bereich Chloridometer ISE-Elektrode | TERMES Instruments Unités Niveau Moyenne Plage attendue Chloromètre Électrode ISE | TERMINI Strumenti Unità Livello Media Range previsto Clorometro Elettrodo ISE | TÉRMINOS Instrumentos Unidades Concentración Media Intervalo esperado Clorímetro Electrodo ISE |
| FOOTNOTES * Mean and range based on limited data N/A Not applicable | FUBNOTEN * Mittelwert und Bereich basieren auf begrenzten Daten N/Z nicht zutreffend | NOTES * Moyennes et plages basées sur des données limitées N/A Non applicable | NOTE A PIÈ DI PAGINA * Media e range basati su dati limitati N/A Non applicabile | NOTAS A PIE DE PÁGINA * Media e intervalo basados en pocos datos. N/P No procede |

European Conformity
CE-Konformitätskennzeichnung
Conformité aux normes européennes
Conformità europea
Conformidad europea

Lot Number
Bezeichnung
Designation du lot
Numero di lotto
Denominación de lote

Manufactured by
Hergestellt von
Fabriqué par
Fabricato da
Fabricado por

For in vitro diagnostic use
In-vitro Diagnosticum
Pour diagnostic in vitro
Per uso diagnostico in vitro
De uso diagnóstico in vitro

Biological Risk
Biogefährdung
Risque biologique
Rischio biologico
Peligro biológico

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Sweat Control for Cystic Fibrosis Testing

Level 1, 2 & 3

ENGLISH

Intended Use - The Quantimetrix Sweat Controls for Cystic Fibrosis Testing are assayed controls intended as a means of monitoring the measurement of conductivity and osmolality by cystic fibrosis screening procedures and to validate the analysis of chloride, potassium, and sodium by confirmatory methods used to quantitate patient samples⁵.

Product Description - The Quantimetrix Sweat Controls for Cystic Fibrosis Testing are liquid, ready to use, requiring no reconstitution or dilution. They are supplied in three levels, nine vials per box, 3 x 3 mL each level. They are prepared in an aqueous simulated human sweat matrix⁶ to which preservatives have been added.

Warnings and Precautions - The Quantimetrix Sweat Controls for Cystic Fibrosis Testing do not contain human source material. No special precautions are required in handling the product other than those routinely used in the laboratory.

Storage and Stability - The Sweat Controls should be stored refrigerated (2 - 8 C) or at room temperature (20 - 25 C). Do not freeze. When stored unopened at 2 - 8 C or at 20 - 25 C the controls are stable until the expiration date stated on the label. When stored at 2 - 8 C between each use the controls are stable for 3 months after opening. Discard the controls if turbid or any evidence of microbial contamination is present.

Procedure - Remove the controls from the refrigerator and allow them to come to room temperature (20 - 25 C), at least 15 minutes, depending on remaining volume. Invert gently to ensure homogeneity of the contents. Treat the controls as you would a patient sample as noted below.

A. Chloride Analysis - If the patient's sweat sample is collected on gauze or filter paper add at least 100 µL of each of the control levels into separate preweighed vials containing the same gauze or filter paper used to collect the patient sample. Treat the controls in the same manner as a patient sample when eluting the control samples and follow the chloridometer manufacturer's instructions for sample titration. Calculate control results in the same manner as patient results⁴.

If the patient's sweat sample is collected in microbore tubing, add 50 µL of each of the control levels into separate titration vials and treat the controls in the same manner as patient samples, following the chloridometer manufacturer's instructions for sample titration or the electrolyte analyzer manufacturer's instructions for sample analysis. Calculate control results in the same manner as patient results.

B. Potassium and Sodium Analysis - If the patient's sweat sample is collected on gauze or filter paper add 50 µL of each of the control levels into separate preweighed vials containing the same gauze or filter paper used to collect the patient sample. Treat the controls in the same manner as a patient sample when eluting the control samples and follow the flame photometer manufacturer's instructions for sample analysis. Calculate control results in the same manner as patient results⁴.

If the patient's sweat sample is collected in microbore tubing, add 10 µL of each of the control levels into separate vials and treat the controls in the same manner as patient samples, following the flame photometer or the electrolyte analyzer manufacturer's instructions for sample analysis. The control and patient results can be read directly. No calculation is necessary.

C. Conductivity Analysis - Treat the controls in the same manner as a patient sample in accordance with the instrument manufacturer's requirements of the test method.

D. Osmolality Analysis - Treat the controls in the same manner as a patient sample in accordance with the instrument manufacturer's requirements of the test method.

Expected Values - The expected values have been established in the Quantimetrix laboratory and from interlaboratory data using the listed manufacturer's instruments. Use of control materials having known component concentrations are an integral part of diagnostic procedures. Monitoring of control values establishes intralaboratory parameters for accuracy and precision of the test method. Each laboratory should establish its own quality control standards for the test method used. A quality control program, Quantrol, is available from Quantimetrix. Please inquire.

Limitations - Any future changes made by the manufacturer of a test method may give a different value from the one indicated. Limitations of the test method are included in the information provided by the instrument manufacturer. Technical updates can be found on our website.

References - ⁴Schales O, Schales SS. A Simple and Accurate Method for the Determination of Chloride in Biological Fluids. J Biol Chem. 1941;140:879-884. ⁵LeGrys VA: Sweat Testing, Clinical Chemistry News, April 1992, page 7-8. ⁶Emrich HM, Stoll E, Colombo JP, Richterrich R, Rossi E: Sweat Composition in Relation to Rate of Sweating in Patients with Cystic Fibrosis of the Pancreas, Pediatric Research, 2: 464-478 (1968).

⁵National Committee for Clinical Laboratory Standards. Sweat Testing: Sample Collection and Quantitative Analysis; Proposed Guideline. NCCLS Document C34-P,

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Bibliographie - *Schales O, Schales SS. A Simple and Accurate Method for the Determination of Chloride in Biological Fluids. J Biol Chem. 1941;140:879-884. *LeGrys VA: Sweat Testing, Clinical Chemistry News, April 1992, page 7-8. *Emrich HM, Stoll E, Colombo JP, Richterrich R, Rossi E: Sweat Composition in Relation to Rate of Sweating in Patients with Cystic Fibrosis of the Pancreas, Pediatric Research, 2: 464-478 (1968). *National Committee for Clinical Laboratory Standards. Sweat Testing: Sample Collection and Quantitative Analysis; Proposed Guideline. NCCLS Document C34-P, NCCLS Villanova PA 19085, 1993.

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Références - *Schales O, Schales SS. A Simple and Accurate Method for the Determination of Chloride in Biological Fluids. J Biol Chem. 1941;140:879-884. *LeGrys VA: Sweat Testing, Clinical Chemistry News, April 1992, page 7-8. *Emrich HM, Stoll E, Colombo JP, Richterrich R, Rossi E: Sweat Composition in Relation to Rate of Sweating in Patients with Cystic Fibrosis of the Pancreas, Pediatric Research, 2: 464-478 (1968). *National Committee for Clinical Laboratory Standards. Sweat Testing: Sample Collection and Quantitative Analysis; Proposed Guideline. NCCLS Document C34-P, NCCLS Villanova PA 19085, 1993.

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Letteratura di riferimento - ¹Schales O, Schales SS. A Simple and Accurate Method for the Determination of Chloride in Biological Fluids. J Biol Chem. 1941;140:879-884. ²LeGrys VA: Sweat Testing, Clinical Chemistry News, April 1992, page 7-8. ³Emrich HM, Stoll E, Colombo JP, Richterrich R, Rossi E: Sweat Composition in Relation to Rate of Sweating in Patients with Cystic Fibrosis of the Pancreas, Pediatric Research, 2: 464-478 (1968). ⁴National Committee for Clinical Laboratory Standards. Sweat Testing: Sample Collection and Quantitative Analysis; Proposed Guideline. NCCLS Document C34-P, NCCLS Villanova PA 19085, 1993.

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Bibliografía - ¹Schales O, Schales SS. A Simple and Accurate Method for the Determination of Chloride in Biological Fluids. J Biol Chem. 1941;140:879-884. ²LeGrys VA: Sweat Testing, Clinical Chemistry News, April 1992, page 7-8. ³Emrich HM, Stoll E, Colombo JP, Richterrich R, Rossi E: Sweat Composition in Relation to Rate of Sweating in Patients with Cystic Fibrosis of the Pancreas, Pediatric Research, 2: 464-478 (1968). ⁴National Committee for Clinical Laboratory Standards. Sweat Testing: Sample Collection and Quantitative Analysis; Proposed Guideline. NCCLS Document C34-P, NCCLS Villanova PA 19085, 1993.