NOTICE: Insert format has changed. • **HINWEIS:** Das Format der Packungsbeilage wurde geändert. • **AVIS:** Le format de la notice d'utilisation a changé.
AVVISO: il formato dell'insero informativo è stato cambiato. • **AVISO:** El formato del prospecto ha cambiado.

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English

Intended Use

The Quantimetrix Sweat Controls for Cystic Fibrosis Testing are assayed controls intended as a means of monitoring the measurement of conductivity and osmolality by cystic fibrosis screening procedures and to validate the analysis of chloride, potassium, and sodium by confirmatory methods used to quantitate patient samples.

Product Description

The Quantimetrix Sweat Controls for Cystic Fibrosis Testing are liquid, ready to use, requiring no reconstitution or dilution. They are supplied in three levels, nine vials per box, 3 x 3 mL each level. They are prepared in an aqueous simulated human sweat matrix to which preservatives have been added.

Warnings and Precautions

The Quantimetrix Sweat Controls for Cystic Fibrosis Testing do not contain human source material. No special precautions are required in handling the product other than those routinely used in the laboratory.

Storage and Stability

The Sweat Controls should be stored refrigerated (2° – 8° C) or at room temperature (20° – 25° C). Do not freeze. When stored unopened at 2° – 8° C or at 20° – 25° C the controls are stable until the expiration date stated on the label. When stored at 2° – 8° C between each use the controls are stable for 3 months after opening. Discard the controls if turbid or any evidence of microbial contamination is present.

Procedure

Remove the controls from the refrigerator and allow them to come to room temperature (20° – 25° C), at least 15 minutes, depending on remaining volume. Invert gently to ensure homogeneity of the contents. Immediately recap the control when not in use.

Expected Values

The expected values have been established in the Quantimetrix laboratory and from interlaboratory data using the listed manufacturer's instruments. Use of control materials having known component concentrations are an integral part of diagnostic procedures. Monitoring of control values establishes intralaboratory parameters for accuracy and precision of the test method. Each laboratory should establish its own quality control standards for the test method used. A quality control program, Quantrol, is available from Quantimetrix. Please inquire.

Limitations

Any future changes made by the manufacturer of a test method may give a different value from the one indicated. Limitations of the test method are included in the information provided by the instrument manufacturer. Technical updates can be found on our website.

References

*Schales O, Schales SS. A Simple and Accurate Method for the Determination of Chloride in Biological Fluids. J Biol Chem. 1941;140:879-884. *LeGrys VA: Sweat Testing, Clinical Chemistry News, April 1992, page 7-8. *Emrich HM, Stoll E, Colombo JP, Richterrich R, Rossi E: Sweat Composition in Relation to Rate of Sweating in Patients with Cystic Fibrosis of the Pancreas, Pediatric Research, 2: 464-478 (1968). *National Committee for Clinical Laboratory Standards. Sweat Testing: Sample Collection and Quantitative Analysis; Proposed Guideline. NCCLS Document C34-A3, NCCLS Wayne PA 19087, 2009.

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¹Schales O, Schales SS. A Simple and Accurate Method for the Determination of Chloride in Biological Fluids. J Biol Chem. 1941;140:879-884. ²LeGrys VA: Sweat Testing, Clinical Chemistry News, April 1992, page 7-8. ³Emrich HM, Stoll E, Colombo JP, Richterrich R, Rossi E: Sweat Composition in Relation to Rate of Sweating in Patients with Cystic Fibrosis of the Pancreas, Pediatric Research, 2: 464-478 (1968). ⁴National Committee for Clinical Laboratory Standards. Sweat Testing: Sample Collection and Quantitative Analysis; Proposed Guideline. NCCLS Document C34-A3, NCCLS Wayne PA 19087, 2009.

Sweat Control for Cystic Fibrosis Testing / Level 1, Level 2, Level 3

| Instruments / Analytes | Units | Level 1 / Lot 19271 | | Level 2 / Lot 19272 | | Level 3 / Lot 19273 | |
|----------------------------------|--------|---------------------|----------------|---------------------|----------------|---------------------|----------------|
| | | mean | expected range | mean | expected range | mean | expected range |
| Advia* | | | | | | | |
| Chloride | mmol/L | 24 | 19 – 29 | 53 | 42 – 63 | 102 | 82 – 123 |
| Potassium | mmol/L | 7.4 | 5.9 – 8.9 | 17.1 | 13.7 – 20.5 | 34.6 | 27.7 – 41.5 |
| Sodium | mmol/L | 25 | 20 – 30 | 54 | 43 – 65 | 107 | 86 – 129 |
| Dimension | | | | | | | |
| Chloride | mmol/L | 29 | 23 – 35 | 58 | 46 – 70 | 111 | 89 – 133 |
| Potassium | mmol/L | 7.3 | 5.8 – 8.8 | 16.7 | 13.4 – 20.1 | 33.6 | 26.9 – 40.3 |
| Sodium | mmol/L | 26 | 21 – 31 | 54 | 43 – 65 | 105 | 84 – 126 |
| Hitachi* | | | | | | | |
| Potassium | mmol/L | 7.0 | 5.6 – 8.4 | 17.4 | 13.9 – 20.9 | 39.6 | 31.7 – 47.5 |
| Sodium | mmol/L | 30 | 24 – 36 | 56 | 45 – 67 | 108 | 86 – 129 |
| Chloridometer* | | | | | | | |
| Chloride | mmol/L | 23 | 18 – 28 | 51 | 41 – 62 | 106 | 85 – 127 |
| Sweat-Chek | | | | | | | |
| Conductivity | mmol/L | 34 | 27 – 40 | 69 | 55 – 83 | 123 | 98 – 147 |
| Advanced Osmometers* | | | | | | | |
| Osmolality | mOs/kg | 61 | 49 – 74 | 140 | 112 – 168 | 278 | 222 – 333 |
| Vapor Pressure Osmometer* | | | | | | | |
| Osmolality | mOs/kg | 55 | 44 – 66 | 132 | 105 – 158 | 265 | 212 – 318 |

Footnotes | Fußnoten | Notes | Note a piè di pagina | Notas a pie de página

English

* Mean and range based on limited data

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This page intentionally left blank.