Dropper® Urine Dipstick Control / Level 1 & 2

LOT 44691
44692

REF 1440-02

2014-10

European Conformity
CE-Konformitätskennzeichnung
Conformité aux normes européennes
Conformità europea
Conformidad europea

Lot Number
Bezeichnung
Designation du lot
Numero di lotto
Denominación de lote

Manufactured by
Hergestellt von
Fabrique par
Fabricato da
Fabricado por

For in vitro diagnostic use
In-vitro Diagnosticum
Pour diagnostic in vitro
Per uso diagnostico in vitro
De uso diagnóstico in vitro

Biological Risk
Biotiefährdung
Risque biologique
Rischio biologico
Peligro biológico

Contents of kit
Inhalt der Packung
Contenu du coffret
Contenido della confezione
Contenido del estuche

Catalog No.
Bestellnr.
N° de catalogue
Catalogo n.
N° de catálogo

See Product Insert
Siehe Packungsbeilage
Voir notice d'utilisation
Vedere il foglio illustrativo del prodotto
Consulte el folleto del producto

Authorized Representative
Bevollmächtigter
Représentant agréé
Rappresentante autorizzato
Representante autorizado

Temperature limitation
Temperaturbegrenzungen
Limites de température
limiti di temperatura
límite de temperatura

Use by (fast day of month)
Verwendbar bis (letzter Tag des Monats)
Utilisable jusqu'à (dernier jour du mois indiqué)
Da utilizzare prima del (ultimo giorno del mese)
Estable hasta (ultimo día del mes)

English

Intended Use

The Quantimetrix Dropper Urine Dipstick Control is intended as a control for urinalysis reagent strips, microalbumin, and creatinine by the listed test methods, and as a control for confirmatory tests such as **Acetest® Clinitest®** and **Ictotest®** reagent tablets, and for **hCG** methods.

Product Description

The Quantimetrix Dropper Urine Dipstick Controls are supplied in two levels, 2 x 25 mL each level per box. They are liquid, ready-to-use, requiring no reconstitution or dilution. They are prepared from human urine fortified to target levels with compounds that produce the desired reaction when tested by the methods indicated in the **Intended Use** section. Preservatives have been added to inhibit microbial growth.

Warnings and Precautions

POTENTIAL BIOHAZARDOUS MATERIAL. Contains human urine and human chorionic gonadotropin (hCG) from pregnancy urine. The human hCG source material used in the manufacture of this product has been tested and found non-reactive for Hepatitis B Surface Antigen, Hepatitis C and HIV 1&2 antibody when tested by FDA accepted methods. No known test method can assure that a product derived from human material does not contain Hepatitis or HIV virus. It is recommended that such samples be handled according to the Centers for Disease Control's Bio-Safety Level 2 recommendations.

Storage and Stability

The Urine Dipstick Control Kit should be stored at 2-8°C when not in use. **Do not freeze.** When stored at 2-8°C the opened Urine Dipstick Control bottles are stable until the expiration date stated on the label. Discard the controls if turbid or any evidence of microbial contamination is present. Discard controls in the same manner as other biological specimens, according to local guidelines.

Procedure for Dipstick Testing

Remove the controls from the refrigerator and allow to come to room temperature (18-25°C), at least 15 minutes, depending on remaining volume. Mix gently by inversion to assure homogeneity of the contents. Avoid foaming. Remove cap and invert bottle. While holding dipstick, gently squeeze the sides of the dropper bottle, and touch the tip of the bottle to the dipstick. Draw across all of the reagent pads, thoroughly saturating each pad. Do not aspirate excess control back into the bottle. Turn dipstick on its side and drain excess control onto absorbent material. Read the urine dipsticks, visually or with an instrumental reader, in accordance with the manufacturers' instructions. Wipe off dropper tips and recap controls. Return them to 2-8°C when not in use.

Procedure for hCG Testing

Note: The bottles of Level 1 Control are to be used as a negative control for hCG methods. The bottles of Level 2 Control are to be used as a positive control for hCG methods.

Note: Most manufacturers of pregnancy test kits specify the volume of sample to be used with their kits. Many kits include transfer pipets to be used to deliver a certain sample volume onto the test device. It is important that sufficient volume be used to produce the correct test result.

The Dropper Control is supplied in dropper tip bottles for convenience and ease of dispensing. Each user should validate that the volume (number of drops) dispensed by the bottle dropper tip is sufficient to meet the pregnancy test kit's requirement for sample volume. Remove the controls from the refrigerator and allow to come to room temperature (18-25°C), at least 15 minutes, depending on remaining volume. Use the hCG positive and negative controls as if they were patient specimens in accordance with the hCG test kit manufacturer's instructions. Immediately recap the controls and return them to 2-8°C when not in use.

Expected Values

For **visual readings**, the expected ranges have been established from interlaboratory data by comparing the dipstick reaction that occurs with the controls to the color comparison chart with multiple lots of each manufacturers' dipsticks or reagent tablets. For expected values for urinalysis reagent strips not listed, please contact Quantimetrix Technical Services. For **instrument readings**, the expected ranges have been established from interlaboratory data from multiple lots of each manufacturers' dipsticks. Each laboratory should establish its own precision parameters. For **specific gravity**, the expected ranges by refractometer have been established from interlaboratory data. For **hCG**, the positive and negative results were obtained by testing each lot number of the controls with multiple lot numbers of different hCG test kits with sensitivities of ≤ 25 mIU/mL.

Limitations

Any future changes made by the manufacturer of a test method may give different values from the indicated range. Detailed information on the limitations of each test method is included in the limitations section of the manufacturers' package insert. Technical updates can be found on our website. The Quality Control Log will no longer be shipped with the product. You can obtain it from the Quantimetrix website at www.quantimetrix.com or contact Tech Support at (310) 536.0006, option 3.

Chemstrip/Combur/DiaScreen/Multistix/Urocheck Users

Colors produced by the **urobilinogen** and/or **bilirubin** reactions on these dipsticks with the Urine Dipstick Control may not be characteristic of those shown on the manufacturer's label when reading the dipstick reactions visually. The urobilinogen reactions are consistent and intensify with the

increase in the urobilinogen concentration but may not provide an exact color match to those displayed on the label.

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Analytes	Level 1 - 44691	Level 2 - 44692
Analyticon® (Visual/Analyzers)		
Ascorbic Acid	Negative - 1+	Negative - 1+
Urobilinogen	Normal	2 - 12 mg/dL (1+ - 4+) ^{9*}
Glucose	Normal	100 - ≥1000 mg/dL (2+ - 5+)
Ketones	Negative	10 - 300 mg/dL (Tr - 3+)
Bilirubin	Negative	1 - 4 mg/dL (1+ - 3+)
Protein	Negative	30 (Tr) - 500 mg/dL
Nitrite	Negative	Positive
Leukocytes	Negative	25 - 500 Leu/μL (1+ - 3+)
Blood	Negative	5 - 300 Ery/μL (1+ - 3+)
pH	5 - 6	6 - 9
Specific Gravity	1.010 - 1.030	1.000 - 1.020
Confirmatory and Other Tests		
Acetest (Ketones)	Negative	Small - Large
Clinitest (Glucose) 5-DROP METHOD	Negative	¼ - ¾
Clinitest (Glucose) 2-DROP METHOD	Negative	½ - 1
Icotest (Bilirubin)	Negative	Positive
Refractometer (Specific Gravity)	1.017 - 1.023	1.008 - 1.016
hCG	Negative	Positive
pH Paper	4 - 6	7 - 9
Sulfosalicylic Acid (Total Protein)	Negative (≤ 0.05 ¹¹)	Positive (≥ 0.50 ¹¹)
DiaScreen® (Visual / Analyzers)		
Urobilinogen	Normal	Not Compatible
Glucose	Negative	50 - 1000 mg/dL ⁸
Ketones	Negative	5 - 160 mg/dL (Tr - 4+)
Bilirubin	Negative	Sm - Lg (1+ - 3+) ^{8,9*}
Protein	Negative	Tr - 300 mg/dL (Tr - 3+)
Nitrite	Negative	Positive
Leukocytes	Negative	Tr - Large (Tr - 3+)
Blood	Negative	Tr - Large (Tr - 3+)
pH	5 - 6	7 - 9
Specific Gravity	1.020 - ≥1.030	1.005 - 1.020
DFI CYBOW/ComboStik/DUS Reagent Strips (Visual)		
Urobilinogen	Normal (0.1~1mg/dL)	2~8mg/dL (33~131umol/L)
Glucose	Negative	50~2000 mg/dL (2.8~111 mmol/L)
Bilirubin	Negative	Small~Large (1+~3+)
Ketones	Negative	5~40mg/dL (0.5~4 mmol/L) (Trace~2+)
Specific Gravity	1.015~1.025	1.005~1.020
Blood	Negative	10~250 RBC/uL (1+~3+)
pH	5~6.5	7~9
Protein, Total	Negative	15~300mg/dL (0.15~3.0g/L) (Trace~3+)
Nitrite	Negative	Positive
Leukocytes	Negative	15~500 WBC/uL (Trace~3+)
Creatinine	10~50 mg/dL (0.9~0.5 mmol/L)	100~300mg/dL (8.8~26.5mmol/L)
Microalbumin	10mg/L	30~150mg/L
DFI CYBOW R-50/ComboStik R-50/DUS R-50		
Bilirubin	Negative	Small~Large (1+~3+)
Urobilinogen	Normal (0.1~1 mg/dL)	2~8 mg/dL (33~131 umol/L)
Ketones	Negative	5~40 mg/dL (0.5~4 mmol/L) (Trace~2+)
Glucose	Negative	50~2000 mg/dL (2.8~111 mmol/L)
Protein	Negative	15~300 mg/dL (0.15~3.0 g/L) (Trace~3+)
Blood	Negative	10~250 RBC/uL (1+~3+)
pH	5~6.5	7~9
Nitrite	Negative	Positive
Leukocytes	Negative	15~500 WBC/uL (Trace~3+)
Specific Gravity	1.015~1.025	1.005~1.020
Micoalbumin	10 mg/L	30~150 mg/L
Creatinine	10~50 mg/dL (0.9~0.5 mmol/L)	100~300 mg/dL (8.8~26.5 mmol/L)
DFI CYBOW Reader 300/ComboStik R-300/DUS R-300/CYBOW Reader 720/ComboStik R-700/DUS R-720		
Urobilinogen	Normal (0.1~1mg/dL)	2~8mg/dL (33~131umol/L)
Glucose	Negative	50~2000mg/dL (2.8~111mmol/L)
Bilirubin	Negative	Small~Large (1+~3+)
Ketones	Negative	5~40mg/dL (0.5~4mmol/L) (Trace~2+)
Specific Gravity	1.015~1.025	1.005~1.020
Blood	Negative	10~250RBC/uL (1+~3+)
pH	5~6.5	7~9
Protein, Total	Negative	15~300mg/dL (0.15~3.0g/L) (Trace~3+)
Nitrite	Negative	Positive
Leukocytes	Negative	15~500 WBC/uL (Trace~3+)
Creatinine	10~50mg/dL (0.9~0.5 mmol/L)	100~300mg/dL (8.8~26.5 mmol/L)
Microalbumin	10mg/L	30~150mg/L

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Analytes	Level 1 - 44691	Level 2 - 44692
Henry Schein® OneStepPro Analyzer* 12		
Blood	Negative	Tr - Lg (Tr - 3+)
Urobilinogen	Normal (0.2 E.U./dL)	1 - ≥8 E.U./dL
Bilirubin	Negative	Sm - Lg (1+ - 3+)
Protein	Negative	Tr - ≥300 mg/dL (Tr - 3+)
Nitrite	Negative	Positive
Ketone	Negative	Tr - ≥160 mg/dL (Tr - 4+)
Glucose	Negative	100 - ≥1000 mg/dL (Tr - 4+)
pH	5 - 6	7 - 9
Specific Gravity	1.010 - 1.030	1.005 - 1.030
Leukocyte	Negative	Tr - Lg (Tr - 3+)
IRIS Diagnostics® (iChem®VELOCITY™ Analyzer)		
Bilirubin	Negative	Not compatible
Urobilinogen	Negative - Normal	2 - ≥8 mg/dL (1+ - 3+)
Ketones	Negative	10 - 150 mg/dL (1+ - 4+)
Glucose	Negative	50 - ≥1000 mg/dL (Tr - 4+)
Protein	Negative	50 - 300 mg/dL (1+ - 3+)
Blood	Negative	0.2 - ≥1 mg/dL (2+ - 3+)
pH	5 - 6	7 - 9
Nitrite	Negative	Positive
Leukocytes	Negative	75 - 500 mg/dL (1+ - 3+)
Specific Gravity	1.000 - 1.020	1.000 - 1.015
MACHERY-NAGEL® (Visual)		
Glucose	Negative - Normal (20 mg/dL)	50 - ≥1000 mg/dL (1+ - 4+) ⁹
Bilirubin	Negative	1 - 4 mg/dL (1+ - 3+) ⁸
Ketones	Negative	25 - 300 mg/dL (1+ - 3+)
Specific Gravity	1.010 - 1.030	1.000 - 1.015
Blood	Negative	5 - 250 Ery/μL (1+ - 3+)
pH	5 - 6	7 - 9
Protein	Negative	30 - 500 mg/dL (1+ - 3+)
Urobilinogen	Normal	2 - 12 mg/dL (1+ - 4+) ⁸
Nitrite	Negative	Positive
Leukocytes	Negative	25 - 500 Leuco/μL (1+ - 3+)
MACHERY-NAGEL® URYXXON® Relax/300/500 Analyzer		
Glucose	Negative - Normal (20 mg/dL)	50 - > 500 mg/dL (1+ - 3+)
Bilirubin	Negative	1 - 4 mg/dL (1+ - 3+) ⁸
Ketones	Negative	25 - 300 mg/dL (1+ - 3+)
Specific Gravity	1.010 - 1.030	1.000 - 1.015
Blood	Negative	10 - 250 Ery/μL (1+ - 3+)
pH	5 - 6	7 - 9
Protein	Negative	30 - 500 mg/dL (1+ - 3+)
Urobilinogen	Normal	2 - 12 mg/dL (1+ - 4+) ⁸
Nitrite	Negative	Positive
Leukocytes	Negative	25 - 500 Leuco/μL (1+ - 3+)
McKesson® (Visual) (10SG Urine Reagent Strips: MFR#121-10SG)¹		
Blood	Negative	Trace - 200 Ery/μL (± - 3+)
Urobilinogen	Normal (0.2 mg/dL)	1 - 12 mg/dL
Bilirubin	Negative	1 - 4 mg/dL (1+ - 3+)
Protein	Negative	15 - 300 mg/dL (± - 3+)
Nitrite	Negative	Positive
Ketone	Negative	5 - 160 mg/dL (± - 4+)
Glucose	Negative	100 - 1000 mg/dL (± - 3+)
pH	5 - 6	7 - 9
Specific Gravity	1.010 - 1.030	1.005 - 1.030
Leukocyte	Negative	15 - 500 Leu/μL (± - 3+)
McKesson® 120 Urine Analyzer ¹		
Blood	Negative	Trace - 200 Ery/μL (± - 3+)
Urobilinogen	Normal (0.2 mg/dL)	Not Compatible
Bilirubin	Negative	1 - 4 mg/dL (1+ - 3+)
Protein	Negative	15 - 300 mg/dL (± - 3+)
Nitrite	Negative	Positive
Ketone	Negative	5 - 160 mg/dL (± - 4+)
Glucose	Negative	100 - 1000 mg/dL (± - 3+)
pH	5 - 6	7 - 9
Specific Gravity	1.015 - 1.030	1.005 - 1.020
Leukocyte	Negative	15 - 500 Leu/μL (± - 3+)

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This page intentionally left blank.