**NOTICE:** Insert format has changed. • **HINWEIS:** Das Format der Packungsbeilage wurde geändert. • **AVIS:** Le format de la notice d'utilisation a changé.  
**AVVISO:** il formato dell'inserito informativo è stato cambiato. • **AVISO:** El formato del prospecto ha cambiado.

# QuantScopics Urine Microscopics Control Level / Level 1 & 2

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## English

### Intended use

The QuantScopics Urine Microscopics Control is intended for use as quality control material for the microscopic evaluation of various urine sediment methods.

Control materials having known component concentrations are an integral part of diagnostic procedures. Daily monitoring of control values establishes intralaboratory parameters for accuracy and precision of the test method. Microscopic QC controls should be run each day the test is performed.<sup>1</sup>

### Summary and Explanation

Standardized microscopic evaluation of urine sediment is an important part of routine analysis of urine. Along with physical and chemical analysis, microscopic examination of urine can provide valuable information regarding not only renal and urinary tract disease, but also metabolic diseases unrelated to the kidney.<sup>2</sup> Urinary sediment microscopy generally includes the detection and identification of red blood cells, leukocytes, epithelial cells, bacteria, casts and crystals.<sup>3,4</sup>

### Product Description

The QuantScopics Urine Microscopics Control is supplied in two levels, 4 x 120 mL bottles per level per box. It is liquid, ready-to-use, requiring no reconstitution or dilution. The control is made from human urine to which stabilized human red blood cells, human white blood cells, and uric acid, calcium oxalate and calcium phosphate crystals have been added. Preservatives have been added to maintain product integrity.

### Warnings and Precautions

**POTENTIAL BIOHAZARDOUS MATERIAL.** All blood donor units comprising the human source material used in the manufacture of this product have been tested and found non-reactive for Hepatitis B surface antigen, Hepatitis C, and HIV antibody when tested by FDA accepted methods. No known test method can assure that a product derived from human blood does not contain Hepatitis or HIV virus. It is recommended that such samples be handled according to the Centers for Disease Control's Bio-Safety Level 2 recommendations.

### Storage and Stability

The QuantScopics Urine Microscopics Control should be stored tightly capped at 2 – 8°C when not in use. Do not freeze. When stored at 2 – 8°C the controls are stable until the expiration date stated on the label. Discard the controls if turbid or any evidence of microbial contamination is present.

### Procedure

Remove the control from the refrigerator and replace the cap on the control bottle with the pour seal cap included in the control box. Allow the control to come to room temperature (18 – 25°C) for approximately 15 – 90 minutes depending on the volume remaining in the bottle. Mix the control thoroughly by inverting the bottle at least 20 times to assure homogeneity of the contents. Avoid foaming. Thorough mixing with each use is important in order to obtain reproducible results. Treat the control as you would a patient sample in accordance with the manufacturer's instructions for the standardized microscopic urinalysis system you are using. The National Committee for Clinical Laboratory Standards (NCCLS) recommends the use of standardized systems in order to yield standardized, reproducible results and to enable the reporting of abnormal sediment elements per unit volume.<sup>1</sup> Immediately close the pour seal cap and store the control at 2 – 8°C when not in use.

### Expected Ranges

The expected ranges for each type of formed element found in this lot of QuantScopics Control is listed in the table on page 4. These ranges were determined by assay of multiple bottles of the indicated lot by the methods listed. A 12 mL sample volume was used for each testing method. The samples were centrifuged at 400 RCF (relative centrifugal force) for 5 minutes. After centrifugation, urine sediment was resuspended in either ~0.5 or ~1.0 mL of remaining supernatant according to the respective cell count method procedure. The ranges listed are based on the range of elements observed in 10 high power fields. Use of other systems or protocols may yield differing results. Each laboratory should establish its own precision parameters.

### Limitations

Any future changes made by the manufacturer of a test system may give different values from the indicated range. Detailed information on the limitations of each test system is included in the limitations section of the system manufacturer's package insert. Technical updates can be found on our website. The Quality Control Log will no longer be shipped with the product. You can obtain it from the Quantimetrix website at [www.4qc.com](http://www.4qc.com) or contact Tech Support at (310) 536-0006, option 3.

### References

<sup>1</sup>NCCLS Document GP16-A. Routine urinalysis and collection, transportation, and preservation of urine specimens; approved guideline. Wayne PA: National Committee for Clinical Laboratory Standards, 1995. <sup>2</sup>Bradley GM, Benson ES. Todd-Sanford Clinical Diagnosis by Laboratory Methods, 15th ed. Philadelphia PA: Saunders, 1974. <sup>3</sup>Bologna CV. Understanding Laboratory Medicine. St. Louis MO: CV Mosby Co, 1971. <sup>4</sup>Dudas HC. Lab Med 1981;12:765.

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# Quantscopics Urine Microscopics Control Level / Level 1 & 2

Microscopic  
Evaluation of  
Urine Sediment

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## Level 1 48241

## Level 2 48242

Analytes / Method	Level 1 48241	Level 2 48242
<b>Red Blood Cells</b>		
COBIO XS, sediMAX, UriSed	0 – 20	20 – 100
IRIS® Diagnostics iQ200 Series Analyzers	3 – 29	17 – 121
Kova® Glasstic® Slide 10 with Grids	3 – 13	29 – 95
Non-grid slides (~0.5 mL resuspended sediment)	0 – 7	2 – 46
Non-grid slides (~1.0 mL resuspended sediment)	0 – 7	2 – 32
Cen-Slide® System	0 – 11	0 – 58
Slide & Coverslip (~0.5 mL resuspended sediment)	0 – 8	2 – 35
Slide & Coverslip (~1.0 mL resuspended sediment)	0 – 7	1 – 20
<b>White Blood Cells</b>		
COBIO XS, sediMAX, UriSed	0 – 10	5 – 50
IRIS Diagnostics iQ200 Series Analyzers	1 – 12	10 – 48
Kova Glasstic Slide 10 with Grids	2 – 13	13 – 54
Non-grid slides (~0.5 mL resuspended sediment)	0 – 8	1 – 40
Non-grid slides (~1.0 mL resuspended sediment)	0 – 5	1 – 18
Cen-Slide System	0 – 11	0 – 36
Slide & Coverslip (~0.5 mL resuspended sediment)	0 – 6	1 – 19
Slide & Coverslip (~1.0 mL resuspended sediment)	0 – 7	1 – 14
<b>Casts</b>		
COBIO XS, sediMAX, UriSed	none	none
IRIS Diagnostics iQ200 Series Analyzers	none	none
Kova Glasstic Slide 10 with Grids	none	none
Non-grid slides (~0.5 mL resuspended sediment)	none	none
Non-grid slides (~1.0 mL resuspended sediment)	none	none
Cen-Slide System	none	none
Slide & Coverslip (~0.5 mL resuspended sediment)	none	none
Slide & Coverslip (~1.0 mL resuspended sediment)	none	none
<b>Crystals</b>		
COBIO XS, sediMAX, UriSed	may be present	present
IRIS Diagnostics iQ200 Series Analyzers	may be present	may be present
Kova Glasstic Slide 10 with Grids	may be present	may be present
Non-grid slides (~0.5 mL resuspended sediment)	may be present	may be present
Non-grid slides (~1.0 mL resuspended sediment)	may be present	may be present
Cen-Slide System	may be present	may be present
Slide & Coverslip (~0.5 mL resuspended sediment)	may be present	may be present
Slide & Coverslip (~1.0 mL resuspended sediment)	may be present	may be present
<b>Bacteria</b>		
COBIO XS, sediMAX, UriSed	may be present	may be present
IRIS Diagnostics iQ200 Series Analyzers	may be present	may be present
Kova Glasstic Slide 10 with Grids	may be present	may be present
Non-grid slides (~0.5 mL resuspended sediment)	may be present	may be present
Non-grid slides (~1.0 mL resuspended sediment)	may be present	may be present
Cen-Slide System	may be present	may be present
Slide & Coverslip (~0.5 mL resuspended sediment)	may be present	may be present
Slide & Coverslip (~1.0 mL resuspended sediment)	may be present	may be present